



ITC R

IN-THE-CANAL RECHARGEABLE

Livio Edge AI 2400 | 2000 | 1600

Livio AI 2400 | 2000 | 1600

Livio 2400 | 2000 | 1600

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone + Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- ITC R 40
- ITC R 50
- ITC R 60
- ITC R 70



Patient Features

- Tinnitus Technology
- Wireless Connectivity
- Rechargeable

Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

ITC R Gain Data

Measurement

	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	110-130	122-139
HFA OSPL90 (dB SPL)	106-124	N/A
RTF OSPL90 (dB SPL)	N/A	106-130
Peak Gain (dB)	40-70	54-79
HFA Full-On Gain (dB)	36-64	N/A
RTF Full-On Gain (dB)	N/A	37-68

Frequency Range (Hz)	<100-9600	<100-9600
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A
Reference Test Gain (dB)	29-47	30-55
Equivalent Input Noise (dB)	<25	<25

Harmonic Distortion

500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3





Estimated Lithium ION Battery Life

Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
-----------------------------------	--------------------	--------------------

Tinnitus Therapy Stimulus

Max RMS Output (dB SPL)	87
Weighted RMS Output Level (dB SPL)	87
Max 1/3 Octave Output (dB SPL)	87

▶ Matrices: 110/40, 115/50, 120/60, 130/70





ITC R

IN-THE-CANAL RECHARGEABLE

Livio Edge AI 1200 | 1000

Livio AI 1200 | 1000

Livio 1200 | 1000

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone +
- Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- ITC R 40
- ITC R 50
- ITC R 60
- ITC R 70



Patient Features

- Tinnitus Technology
- Wireless Connectivity
- Rechargeable

Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

ITC R Gain Data

Measurement

	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	110-130	122-139
HFA OSPL90 (dB SPL)	106-124	N/A
RTF OSPL90 (dB SPL)	N/A	106-130
Peak Gain (dB)	40-70	54-79
HFA Full-On Gain (dB)	36-64	N/A
RTF Full-On Gain (dB)	N/A	37-68

Frequency Range (Hz)	<100-7900	<100-8000
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A
Reference Test Gain (dB)	29-47	30-55
Equivalent Input Noise (dB)	<25	<25

Harmonic Distortion

500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3





Estimated Lithium ION Battery Life

Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
-----------------------------------	--------------------	--------------------

Tinnitus Therapy Stimulus

Max RMS Output (dB SPL)	87
Weighted RMS Output Level (dB SPL)	87
Max 1/3 Octave Output (dB SPL)	87

▶ Matrices: 110/40, 115/50, 120/60, 130/70**



*Results will vary based on wireless usage.

**The highest matrix of 130/70 is not offered on 1000 tier products.



HS R

HALF SHELL RECHARGEABLE

Livio Edge AI 2400 | 2000 | 1600

Livio AI 2400 | 2000 | 1600

Livio 2400 | 2000 | 1600

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone +
- Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- HS R 40
- HS R 50
- HS R 60
- HS R 70



Patient Features

- Tinnitus Technology
- Wireless Connectivity
- Rechargeable

Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

HS R Gain Data

Measurement

	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	110-130	122-139
HFA OSPL90 (dB SPL)	106-124	N/A
RTF OSPL90 (dB SPL)	N/A	106-130
Peak Gain (dB)	40-70	54-79
HFA Full-On Gain (dB)	36-64	N/A
RTF Full-On Gain (dB)	N/A	37-68

Frequency Range (Hz)	<100-9600	<100-9600
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A
Reference Test Gain (dB)	29-47	30-55
Equivalent Input Noise (dB)	<25	<25

Harmonic Distortion

500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3





Estimated Lithium ION Battery Life

Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
-----------------------------------	--------------------	--------------------

Tinnitus Therapy Stimulus

Max RMS Output (dB SPL)	87
Weighted RMS Output Level (dB SPL)	87
Max 1/3 Octave Output (dB SPL)	87

▶ Matrices: 110/40, 115/50, 120/60, 130/70





HS R

HALF SHELL RECHARGEABLE

Livio Edge AI 1200 | 1000

Livio AI 1200 | 1000

Livio 1200 | 1000

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone +
- Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- HS R 40
- HS R 50
- HS R 60
- HS R 70



Patient Features

- Tinnitus Technology
- Wireless Connectivity
- Rechargeable

Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

HS R Gain Data

Measurement

	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	110-130	122-139
HFA OSPL90 (dB SPL)	106-124	N/A
RTF OSPL90 (dB SPL)	N/A	106-130
Peak Gain (dB)	40-70	54-79
HFA Full-On Gain (dB)	36-64	N/A
RTF Full-On Gain (dB)	N/A	37-68

Frequency Range (Hz)	<100-7900	<100-8000
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A
Reference Test Gain (dB)	29-47	30-55
Equivalent Input Noise (dB)	<25	<25

Harmonic Distortion

500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3





Estimated Lithium ION Battery Life

Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
-----------------------------------	--------------------	--------------------

Tinnitus Therapy Stimulus

Max RMS Output (dB SPL)	87
Weighted RMS Output Level (dB SPL)	87
Max 1/3 Octave Output (dB SPL)	87

▶ Matrices: 110/40, 115/50, 120/60, 130/70**



*Results will vary based on wireless usage.

**The highest matrix of 130/70 is not offered on 1000 tier products.



ITE R

IN-THE-EAR RECHARGEABLE

Livio Edge AI 2400 | 2000 | 1600

Livio AI 2400 | 2000 | 1600

Livio 2400 | 2000 | 1600

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone +
- Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- ITE R 50
- ITE R 60
- ITE R 70



Patient Features

- Tinnitus Technology
- Wireless Connectivity
- Rechargeable
- Telecoil




Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

ITE R Gain Data

Measurement	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	115-130	124-137
HFA OSPL90 (dB SPL)	109-123	N/A
RTF OSPL90 (dB SPL)	N/A	110-130
Peak Gain (dB)	50-70	60-77
HFA Full-On Gain (dB)	46-64	N/A
RTF Full-On Gain (dB)	N/A	46-68
Frequency Range (Hz)	<100-8200	<100-8900
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0,1.6,2.5	N/A
Reference Test Gain (dB)	32-46	35-55
Equivalent Input Noise (dB)	<25	<25
Harmonic Distortion		
500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3
Induction Coil Sensitivity		
HFA SPLITS (ANSI) (dB SPL)	95-108	N/A
MASL (IEC) (dB SPL)	N/A	83-106
Estimated Lithium ION Battery Life		
Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
Tinnitus Therapy Stimulus		
Max RMS Output (dB SPL)	87	
Weighted RMS Output Level (dB SPL)	87	
Max 1/3 Octave Output (dB SPL)	87	

▶ Matrices: 115/50, 120/60, 130/70





ITE R

IN-THE-EAR RECHARGEABLE

Livio Edge AI 1200 | 1000

Livio AI 1200 | 1000

Livio 1200 | 1000

Color Guide

Faceplate Colors



Shell Colors




Accessory Compatibility

- TV
- Remote Microphone +
- Remote
- Mini Remote Microphone
- Table Microphone
- 2.4 GHz Programmer

Fitting Range

- ITE R 50
- ITE R 60
- ITE R 70






Livio Edge AI/ Livio AI Technology

- Healthable hearing technology with embedded sensors and artificial intelligence

ITE R Gain Data

Measurement	ANSI/IEC 2cc Coupler	IEC OES Coupler
Peak OSPL90 (dB SPL)	115-130	124-137
HFA OSPL90 (dB SPL)	109-123	N/A
RTF OSPL90 (dB SPL)	N/A	110-130
Peak Gain (dB)	50-70	60-77
HFA Full-On Gain (dB)	46-64	N/A
RTF Full-On Gain (dB)	N/A	46-68
Frequency Range (Hz)	<100-7700	<100-7800
Reference Test Freq. (kHz)	N/A	1.6
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A
Reference Test Gain (dB)	32-46	35-55
Equivalent Input Noise (dB)	<25	<25
Harmonic Distortion		
500 Hz (%)	<3	<3
800 Hz (%)	<3	<3
1600 Hz (%)	<3	<3
Induction Coil Sensitivity		
HFA SPLITS (ANSI) (dB SPL)	95-108	N/A
MASL (IEC) (dB SPL)	N/A	83-106
Estimated Lithium ION Battery Life		
Li-Ion Rechargeable Battery (hrs)	Up to 24 hours*	Up to 24 hours*
Tinnitus Therapy Stimulus		
Max RMS Output (dB SPL)	87	
Weighted RMS Output Level (dB SPL)	87	
Max 1/3 Octave Output (dB SPL)	87	

▶ Matrices: 115/50, 120/60, 130/70**



*Results will vary based on wireless usage.

**The highest matrix of 130/70 is not offered on 1000 tier products.